



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

5-HT-5A Antibody [#C12018]

Catalog Number: C12018

Concentration: 1mg/ml

Swiss-Prot No.: P47898

Other Names: 5-hydroxytryptamine receptor 5A, 5-HT-5A, Serotonin receptor 5A, 5-HT-5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human 5-HT-5A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 5-HT-5A antibody detects endogenous levels of total 5-HT-5A protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Rees S., FEBS Lett. 355:242-246(1994).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Sjoebloom T., nce 314:268-274(2006).

