



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACTN1 Antibody [#C12025]

Catalog Number: C12025

Concentration: 1mg/ml

Swiss-Prot No.: P12814

Other Names: actinin 1 smooth muscle; actinin, alpha 1; ACTN1; alpha-actinin 1; Alpha-actinin cytoskeletal isoform; Alpha-actinin-1; F-actin cross-linking protein; FLJ40884; FLJ54432; Non-muscle alpha-actinin-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ACTN1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ACTN1 antibody detects endogenous levels of total ACTN1 protein.

Reactivity: Human; Mouse; Rat


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Millake D.B., Nucleic Acids Res. 17:6725-6725(1989).

Youssoufian H., Am. J. Hum. Genet. 47:62-71(1990).

Thorsen K., Mol. Cell. Proteomics 7:1214-1224(2008).



For Research Use Only