



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADCY4 Antibody [#C12034]

Catalog Number: C12034

Concentration: 1mg/ml

Swiss-Prot No.: Q8NFM4

Other Names: Adenylate cyclase, type IV; Adenylyl cyclase 4; ATP pyrophosphate-lyase 4; EC 4.6.1.1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADCY4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADCY4 antibody detects endogenous levels of total ADCY4 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ludwig M.G., J. Recept. Signal Transduct. Res. 22:79-110(2002).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Ota T., Nat. Genet. 36:40-45(2004).

