



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ADCY5 Antibody [#C12035]

**Catalog Number:** C12035

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95622/O43306

**Other Names:** ATP pyrophosphate-lyase 5, Adenylate cyclase type 5, Adenylate cyclase type V, Adenylyl cyclase 5

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ADCY5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ADCY5 antibody detects endogenous levels of total ADCY5 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Muzny D.M., Nature 440:1194-1198(2006).

Ludwig M.G., J. Recept. Signal Transduct. Res. 22:79-110(2002).

Raimundo S., Clin. Chim. Acta 285:155-161(1999).

