



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADD3 Antibody [#C12042]

Catalog Number: C12042

Concentration: 1mg/ml

Swiss-Prot No.: Q9UEY8

Other Names: ADDG, ADDL, adducin, gamma, adducin-like protein 70, gamma adducin, protein kinase C binding protein 35H

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADD3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADD3 antibody detects endogenous levels of total ADD3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Katagiri T., Cytogenet. Cell Genet. 74:90-95(1996).

Citterio L., Biochem. Biophys. Res. Commun. 266:110-114(1999).

Gocke C.B., J. Biol. Chem. 280:5004-5012(2005).

