



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## AIFM2 Antibody [#C12047]

**Catalog Number:** C12047

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BRQ8

**Other Names:** Apoptosis-inducing factor 2, EC 1.-.-.-, Apoptosis-inducing factor-like mitochondrion-associated inducer of death, Apoptosis-inducing factor homologous mitochondrion-associated inducer of death, p53-responsive gene 3 protein, AIFM2, AMID, PRG3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human AIFM2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** AIFM2 antibody detects endogenous levels of total AIFM2 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Ohiro Y., FEBS Lett. 524:163-171(2002).

Wu M., J. Biol. Chem. 277:25617-25623(2002).

Ota T., Nat. Genet. 36:40-45(2004).

