



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AIM2 Antibody [#C12049]

Catalog Number: C12049

Concentration: 1mg/ml

Swiss-Prot No.: O14862

Other Names: Interferon-inducible protein AIM2, Absent in melanoma 2, AIM2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AIM2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AIM2 antibody detects endogenous levels of total AIM2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

DeYoung K.L., Oncogene 15:453-457(1997).

Gregory S.G., Nature 441:315-321(2006).

