



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

ARRC Antibody [#C12058]

Catalog Number: C12058

Concentration: 1mg/ml

Swiss-Prot No.: P36575

Other Names: Arrestin-C, Cone arrestin, cArr, Retinal cone arrestin-3, X-arrestin, C-arrestin, ARR3, ARRX, CAR

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ARRC.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARRC antibody detects endogenous levels of total ARRC protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Murakami A., FEBS Lett. 334:203-209(1993).

Craft C.M., J. Biol. Chem. 269:4613-4619(1994).

Sakuma H., FEBS Lett. 382:105-110(1996).

