



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AVEN Antibody [#C12070]

Catalog Number: C12070

Concentration: 1mg/ml

Swiss-Prot No.: Q9NQS1

Other Names: PDCD12; programmed cell death 12, cell death regulator Aven

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human AVEN.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AVEN antibody detects endogenous levels of total AVEN protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Chau B.N., Mol. Cell 6:31-40(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

