



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH11 Antibody [#C12087]

Catalog Number: C12087

Concentration: 1mg/ml

Swiss-Prot No.: P55287

Other Names: CAD11; cadherin 11, type 2; cadherin 11, type 2, OB-cadherin (osteoblast); cadherin-11; CDHOB; OB; OB-cadherin; OSF-4; osteoblast-cadherin

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CDH11.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH11 antibody detects endogenous levels of total CDH11 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Tanihara H., Cell Adhes. Commun. 2:15-26(1994).

Okazaki M., J. Biol. Chem. 269:12092-12098(1994).

Suzuki S., Cell Regul. 2:261-270(1991).

