



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CADM2 Antibody [#C12109]

Catalog Number: C12109

Concentration: 1mg/ml

Swiss-Prot No.: Q8N3J6

Other Names: Cell adhesion molecule 2 [Precursor], Immunoglobulin superfamily member 4D, Nectin-like protein 3, CADM2, IGSF4D, NECL3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CADM2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CADM2 antibody detects endogenous levels of total CADM2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Gingrich J.R., D'Angelo A., Submitted (AUG-2002) to the EMBL/GenBank/DDBJ databases.

The German cDNA consortium, Submitted (SEP-2004) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

