



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CADM3 Antibody [#C12110]

Catalog Number: C12110

Concentration: 1mg/ml

Swiss-Prot No.: Q8N126

Other Names: Cell adhesion molecule 3 [Precursor], Immunoglobulin superfamily member 4B, Nectin-like protein 1, TSLC1-like protein 1, Synaptic cell adhesion molecule 3, Brain immunoglobulin receptor, CADM3, IGSF4B, NECL1, SYNCAM3, TSLL1, UNQ225/PRO258

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human CADM3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CADM3 antibody detects endogenous levels of total CADM3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:1000 ELISA: 1:5000

References:

Fukuhara H., *Oncogene* 20:5401-5407(2001).

Zhou Y., *Biochim. Biophys. Acta* 1669:142-154(2005).

Clark H.F., *Genome Res.* 13:2265-2270(2003).

