



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CARD11 Antibody [#C12113]

**Catalog Number:** C12113

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BXL7

**Other Names:** Caspase recruitment domain-containing protein 11, CARD-containing MAGUK protein 3, Carma 1, CARD11, CARMA1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

CARD11.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CARD11 antibody detects endogenous levels of total CARD11 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Bertin J., J. Biol. Chem. 276:11877-11882(2001).

Gaide O., FEBS Lett. 496:121-127(2001).

Hillier L.W., Nature 424:157-164(2003).

