



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CDCA2 Antibody [#C12173]

**Catalog Number:** C12173

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q69YH5

**Other Names:** FLJ25804, MGC129906, MGC129907,  
cell division cycle associated 2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from C-terminal of human  
CDCA2.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** CDCA2 antibody detects endogenous levels  
of total CDCA2 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Walker M.G., Curr. Cancer Drug Targets 1:73-83(2001).

Ota T., Genet. 36:40-45(2004).

Olsen J.V., Cell 127:635-648(2006).

