



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## OR8B4 Antibody [#G929]

**Catalog Number:** G929

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96RC9

**Other Names:** Olfactory receptor 8B4; OR8B4P

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human OR8B4.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR8B4 antibody detects endogenous levels of total OR8B4 protein.


**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:5000


**References:**

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Fuchs T., Genomics 80:295-302(2002).



Western blot analysis of extracts from MCF-7 cells and HeLa cells, using OR8B4 antibody.



Immunofluorescence analysis of MCF-7 cells, using OR8B4 antibody.