



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## U12 Antibody [#G956]

**Catalog Number:** G956

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P52381

**Other Names:** G-protein coupled receptor homolog U12

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human U12.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


**Specificity:** U12 antibody detects endogenous levels of total U12 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Nicholas J., Submitted (JAN-1996) to the EMBL/GenBank/DDBJ databases.



**For Research Use Only**