




Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACRO (heavy chain, Cleaved-Ile43) Antibody [#L0125]

- Catalog Number:** L0125
Amount: 50µg/50µl, 100µg/100µl, 200µg/200µl
Swiss-Prot No. : P10323
All Names: Acrosin, EC=3.4.21.10
All Sites: Human: Ile43
Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage/Stability: Store at -20°C/1 year
Immunogen: The antiserum was produced against synthesized peptide derived from human ACRO.
Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Specificity/Sensitivity: ACRO (heavy chain, Cleaved-Ile43) antibody detects endogenous levels of fragment of activated ACRO resulting from cleavage adjacent to Ile43.
Reactivity: Human
Applications: WB: 1:500~1:3000 ELISA: 1:1000
References: Baba T., FEBS Lett. 244:296-300(1989) [PubMed: 2493394].
Keime S., Eur. J. Biochem. 190:195-200(1990) [PubMed: 2114285].
Adham I.M., Hum. Genet. 84:125-128(1990) [PubMed: 2298447]



For Research Use Only