



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AZI2 Antibody [#C11698]

Catalog Number: C11698

Concentration: 1mg/ml

Swiss-Prot No.: Q9H6S1

Other Names: 5-azacytidine-induced protein 2,
NF-kappa-B-activating kinase-associated protein 1,
Nak-associated protein 1, TILP, AZI2, NAP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from N-terminal of human

AZI2.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: AZI2 antibody detects endogenous levels of
total AZI2 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fujita F., Mol. Cell. Biol. 23:7780-7793(2003).

Ota T., Nat. Genet. 36:40-45(2004).

Sasai M., J. Immunol. 174:27-30(2005).



For Research Use Only