



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## DUSP22 Antibody [#C11753]

**Catalog Number:** C11753

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9NRW4

**Other Names:** Dual specificity protein phosphatase 22, EC 3.1.3.48, EC 3.1.3.16, JNK-stimulatory phosphatase-1, JSP-1, Mitogen-activated protein kinase phosphatase x, LMW-DSP2, DUSP22, JSP1, MKPX

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human DUSP22.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** DUSP22 antibody detects endogenous levels of total DUSP22 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Shen Y., Proc. Natl. Acad. Sci. U.S.A. 98:13613-13618(2001).

Mungall A.J., Nature 425:805-811(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

