



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ASF1B Antibody [#C11765]

Catalog Number: C11765

Concentration: 1mg/ml

Swiss-Prot No.: Q9NVP2

Other Names: Histone chaperone ASF1B, Anti-silencing function protein 1 homolog B, hAsf1, hAsf1b, CCG1-interacting factor A-II, CIA-II, hCIA-II, ASF1B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ASF1B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ASF1B antibody detects endogenous levels of total ASF1B protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Sillje H.H.W., Curr. Biol. 11:1068-1073(2001).

ehara T., J. Biol. Chem. 278:35660-35667(2003).

Ota T., Nat. Genet. 36:40-45(2004).

