



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DOK5 Antibody [#C11799]

Catalog Number: C11799

Concentration: 1mg/ml

Swiss-Prot No.: Q9P104

Other Names: Docking protein 5; Downstream of tyrosine kinase 5; Insulin receptor substrate 6; IRS-6; IRS6; DOK5; C20orf180

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DOK5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DOK5 antibody detects endogenous levels of total DOK5 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Luo W.Q., Submitted (MAR-1999) to the EMBL/GenBank/DDBJ databases.

Favre C., Genes Immun. 4:40-45(2003).

Cai D., J. Biol. Chem. 278:25323-25330(2003).

