



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Collagen XXV Antibody [#C12227]

Catalog Number: C12227

Concentration: 1mg/ml

Swiss-Prot No.: Q9BXS0

Other Names: Collagen alpha-1(XXV) chain, CLAC-P, Alzheimer disease amyloid-associated protein, AMY, Collagen-like Alzheimer amyloid plaque component, COPA1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

Collagen XXV.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen XXV antibody detects endogenous levels of total Collagen XXV protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Hashimoto T., EMBO J. 21:1524-1534(2002).

Soederberg L., J. Neuropathol. Exp. Neurol. 62:1108-1117(2003).

Soederberg L., J. Biol. Chem. 280:1007-1015(2005).

