



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

COX7S/A2 Antibody [#C12243]

Catalog Number: C12243

Concentration: 1mg/ml

Swiss-Prot No.: O60397

Other Names: Cytochrome c oxidase subunit VIIa 3, mitochondrial [Precursor], COX7AP2, COX7A3, COX7AL2, COX7S

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

COX7S/A2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: COX7S/A2 antibody detects endogenous levels of total COX7S/A2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:10000


References:

Hillier L.W., Nature 434:724-731(2005).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using COX7S/A2 antibody.



Western blot analysis of extracts from rat heart cells, using COX7S/A2 antibody.