



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Cytochrome P450 2C8 Antibody [#C12262]

Catalog Number: C12262

Concentration: 1mg/ml

Swiss-Prot No.: P10632

Other Names: Cytochrome P450 2C8, EC 1.14.14.1, CYP11C8, P450 form 1, P450 MP-12/MP-20, P450 IIC2, S-mephenytoin 4-hydroxylase, CYP2C8, CP2C8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

Cytochrome P450 2C8.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 2C8 antibody detects endogenous levels of total Cytochrome P450 2C8 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Romkes M., Biochemistry 30:3247-3255(1991).

Ged C., Biochim. Biophys. Acta 1088:433-435(1991).

Zeldin D.C., Arch. Biochem. Biophys. 322:76-86(1995).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human testis tissue, using Cytochrome P450 2C8 antibody.