



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Cytochrome P450 2D6 Antibody [#C12266]

**Catalog Number:** C12266

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P10635

**Other Names:** Cytochrome P450 2D6, EC 1.14.14.1, CYP11D6, P450-DB1, Debrisoquine 4-hydroxylase, CYP2D6, CYP2DL1, CP2D6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 2D6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Cytochrome P450 2D6 antibody detects endogenous levels of total Cytochrome P450 2D6 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Tyndale R., Pharmacogenetics 1:26-32(1991).

Yokota H., Pharmacogenetics 3:256-263(1993).

Evert B., Naunyn Schmiedebergs Arch. Pharmacol. 350:434-439(1994).

