



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Cytochrome P450 2U1 Antibody [#C12272]

**Catalog Number:** C12272

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q7Z449

**Other Names:** Cytochrome P450 2U1, EC 1.14.14.1, CP2U1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 2U1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Cytochrome P450 2U1 antibody detects endogenous levels of total Cytochrome P450 2U1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Chuang S.S., J. Biol. Chem. 279:6305-6314(2004).

Karlgren M., Biochem. Biophys. Res. Commun. 315:679-685(2004).

Choudhary D., Arch. Biochem. Biophys. 436:50-61(2005).

