



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## BAX Antibody [#C13029]

**Catalog Number:** C13029

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P55269

**Other Names:** BAX protein, cytoplasmic isoform delta, BAX

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human BAX.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** BAX antibody detects endogenous levels of total BAX protein.


**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%);  
Rat (Identities = 100%, Positives = 100%)

**Applications:** IHC: 1:50~1:100 IF: 1:100~1:500  
ELISA: 1:1000

**References:**


Apte S.S., Genomics 26:592-594(1995).

Guo B., Nature 423:456-461(2003).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using BAX antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using BAX antibody.