



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL15 Antibody [#C14059]

Catalog Number: C14059

Concentration: 1mg/ml

Swiss-Prot No.: Q9P015

Other Names: 39S ribosomal protein L15, mitochondrial [Precursor]; L15mt; MRP-L15; HSPC145

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPL15.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL15 antibody detects endogenous levels of total MRPL15 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Zhang Q.-H., Genome Res. 10:1546-1560(2000).

Wiemann S., Genome Res. 11:422-435(2001).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

