



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL18 Antibody [#C14062]

Catalog Number: C14062

Concentration: 1mg/ml

Swiss-Prot No.: Q9H0U6

Other Names: 39S ribosomal protein L18, mitochondrial; L18mt; mitochondrial ribosomal protein L18; MRP-L18; RM18

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRPL18.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL18 antibody detects endogenous levels of total MRPL18 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Wiemann S., Genome Res. 11:422-435(2001).

Zhang Q.-H., Genome Res. 10:1546-1560(2000).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

