



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRPL20 Antibody [#C14065]

Catalog Number: C14065

Concentration: 1mg/ml

Swiss-Prot No.: Q9BYC9

Other Names: 39S ribosomal protein L20, mitochondrial [Precursor]; L20mt; MRP-L20; MRPL20

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human MRPL20.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRPL20 antibody detects endogenous levels of total MRPL20 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Suzuki T., J. Biol. Chem. 276:21724-21736(2001).

Gregory S.G., Nature 441:315-321(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

