



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## RS27L Antibody [#C14115]

**Catalog Number:** C14115

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q71UM5

**Other Names:** 40S ribosomal protein S27-like protein;  
RPS27L

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human RS27L.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** RS27L antibody detects endogenous levels of total RS27L protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Zhang Q.-H., Genome Res. 10:1546-1560(2000).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

