



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABHD12B Antibody [#C14210]

Catalog Number: C14210

Concentration: 1mg/ml

Swiss-Prot No.: Q7Z5M8

Other Names: Abhydrolase domain-containing protein 12B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABHD12B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABHD12B antibody detects endogenous levels of total ABHD12B protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Sjoeblom T., Science 314:268-274(2006).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using ABHD12B antibody.