



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABHD14B Antibody [#C14213]

Catalog Number: C14213

Concentration: 1mg/ml

Swiss-Prot No.: Q96IU4

Other Names: Abhydrolase domain-containing protein 14B; EC 3.-.-.; CCG1-interacting factor B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABHD14B.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: ABHD14B antibody detects endogenous levels of total ABHD14B protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 86%, Positives = 93%);
Rat (Identities = 93%, Positives = 100%)

Applications: IHC: 1:50~1:100 ELISA: 1:5000

References:

Zilfou JT, et al. (2001) Mol Cell Biol 21(12): 3974-3985.
Boehm M, et al. (2002) EMBO J 21(13): 3390-3401.
Vadlamudi RK, et al. (2005) Mol Cell Biol 25(9): 3726-3736.



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ABHD14B antibody.