



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACHE Antibody [#C14230]

Catalog Number: C14230

Concentration: 1mg/ml

Swiss-Prot No.: P22303

Other Names: Acetylcholinesterase; AChE

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ACHE.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: ACHE antibody detects endogenous levels of total ACHE protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Soreq H., Proc. Natl. Acad. Sci. U.S.A. 87:9688-9692(1990).

Karpel R., Exp. Cell Res. 210:268-277(1994).

Yang L., Submitted (JAN-2001) to the EMBL/GenBank/DDBJ databases.

