



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARL2BP Antibody [#C14350]

Catalog Number: C14350

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y2Y0

Other Names: ADP-ribosylation factor-like 2 binding; AR2BP; Arf-like 2 binding protein BART1; ARL2BP; BART1; binder of Arl Two

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ARL2BP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARL2BP antibody detects endogenous levels of total ARL2BP protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Sharer J.D., J. Biol. Chem. 274:27553-27561(1999).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Van Valkenburgh H., J. Biol. Chem. 276:22826-22837(2001).

