



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HIBADH Antibody [#C14092]

Catalog Number: C14092

Concentration: 1mg/ml

Swiss-Prot No.: P31937

Other Names: 3-hydroxyisobutyrate dehydrogenase, mitochondrial; EC 1.1.1.31; HIBADH

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human HIBADH.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HIBADH antibody detects endogenous levels of total HIBADH protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 92%, Positives = 100%);
Rat (Identities = 92%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Hillier L.W., Nature 424:157-164(2003).

Hughes G.J., Electrophoresis 14:1216-1222(1993).

Structural genomics consortium (SGC); Submitted (APR-2006) to the PDB data bank.



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using HIBADH antibody.