



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ATP5L2 Antibody [#C14604]

Catalog Number: C14604

Concentration: 1mg/ml

Swiss-Prot No.: Q7Z4Y8

Other Names: ATP synthase subunit g 2, mitochondrial; ATPase subunit g 2; ATP5K2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ATP5L2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ATP5L2 antibody detects endogenous levels of total ATP5L2 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000


References:

Lin W., Submitted (SEP-1998) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from A549 cells, using ATP5L2 antibody.



Immunofluorescence analysis of A549 cells, using ATP5L2 antibody.