



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCA13 Antibody [#C14611]

Catalog Number: C14611

Concentration: 1mg/ml

Swiss-Prot No.: Q86UQ4

Other Names: ATP-binding cassette sub-family A member 13

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ABCA13.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABCA13 antibody detects endogenous levels of total ABCA13 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Prades C., Cytogenet. Genome Res. 98:160-168(2002).

Ota T., Nat. Genet. 36:40-45(2004).

Schaap F.G., Submitted (APR-2002) to the EMBL/GenBank/DDBJ databases.

