



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ABCB10 Antibody [#C14619]

Catalog Number: C14619

Concentration: 1mg/ml

Swiss-Prot No.: Q9NRK6

Other Names: ABC transporter 10; ABCBA; ATP-binding cassette sub-family B member 10, mitochondrial; ATP-binding cassette transporter 10; M-ABC2; mitochondrial ATP-binding cassette 2; MTABC2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ABCB10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ABCB10 antibody detects endogenous levels of total ABCB10 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Zhang F., FEBS Lett. 478:89-94(2000).

Ito K., Submitted (APR-1998) to the EMBL/GenBank/DDBJ databases.

Gregory S.G., Nature 441:315-321(2006).

