



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ABCF1 Antibody [#C14629]

**Catalog Number:** C14629

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q8NE71

**Other Names:** ABC27; ABC50; ABCF1; ATP-binding cassette 50; ATP-binding cassette 50 (TNF-alpha stimulated); ATP-binding cassette sub-family F member 1; ATP-binding cassette, sub-family F (GCN20), member 1; TNF-alpha-stimulated ABC protein; TNFalpha-inducible ATP-binding protein

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human ABCF1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ABCF1 antibody detects endogenous levels of total ABCF1 protein.

**Reactivity:** Human; Mouse; Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Richard M., Genomics 53:137-145(1998).

Shiina S., Submitted (SEP-1999) to the EMBL/GenBank/DDBJ databases.

Shiina T., Submitted (JUL-2002) to the EMBL/GenBank/DDBJ databases.

